# ALSATOM SU-MPC

ELETTROBISTURI PER CHIRURGIA MONOPOLARE E BIPOLEARE  
DIATHERMY UNITS FOR MONOPOLAR AND BIPOLAR SURGERY

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The new ALSATOM SU-MPC series is composed of intuitive Electrosurgical Units, with very high performances. They have both standard and pulsed currents, which can minimize the thermal effects into the tissues and reduce the harmful smokes, which are produced by all the Electrosurgical Units. By setting these currents, the operators can perform very fine cuts, very similar to the ones obtained with Radio Frequency units, as well as very delicate coagulations, which cannot be reached in any other way. Moreover, these new Electrosurgical Units are provided with a current for Micro-Coagulation, which allows delivering single pulses from 0.1sec to 1sec.

- There are 5 available models:
  - 4 for the Monopolar and Bipolar use, and the Monopolar use under liquid with mini-resectors and needles of 5Fr.
  - 1 for the Bipolar use only with cut and coagulation, as well as cut, coagulation and vaporization in Saline Solution with mini-resectors, needles of 5Fr and Arthroscopy instruments.
- They are controlled through a microprocessor with self-check at the activation, interruption in case of failures "Error codes", and automatic self-regulation of the currents according to the tissues. They do not need any periodic calibration.
- The front panel has got buttons with digital displays. The settings of the cut and coagulation currents are independent (micro-adjustment of the power: 1W from 0 to 30W, and then 2W), and they are stored and re-established at the moment of the switching-on.
- The safety circuit of the neutral electrode also controls the good contact on the tissues, by connecting the ones with the double section.
- The footswitch pedals do not have any electric current, and they work with the air. They can be used into the operating theatre; they are waterproof and anti-explosion. The models SU 140/D-MPC and SU 140/BD-MPC have got the very practical double pedal, in order to activate the cut and coagulation currents.
- The Monopolar currents for the cut and coagulation can also be activated, into the models SU 140-MPC and SU 140/D-MPC, through the handles with two buttons. This mode allows the free activation of the unit: the Monopolar currents by the handle, and the Bipolar coagulation by the footswitch.
- The units can be used in Gynaecology, Dermatology, Plastic and Aesthetic Surgery, ENT, Maxillo Facial Surgery, Dental Surgery, Ophthalmology, Gastroenterology and Veterinary (even in Laparoscopy) and in all the other surgeries (i.e. in Urologic Coagulation under liquid or in Emergency).

## THE CURRENTS OF THE MODELS FOR MONOPOLAR AND BIPOLAR USE

### MONOPOLAR CURRENTS

<b>PURE</b>	Pure Cut, suitable also for the use under liquid in small Hysteroscopy.
<b>P PULSED</b>	Pure Pulsed Cut, suitable for very fine cuts, with a minimum thermal effect (i.e. for conizations in Gynaecology or Blepharoplasty); it is used also to reduce the smokes.
<b>BLEND</b>	Coagulating Cut.
<b>B PULSED 1</b>	Coagulating Pulsed Cut. Similar to BLEND, but indicated to reduce the thermal effect and the smokes.
<b>B PULSED 2</b>	Slow Coagulating Pulsed Cut. Similar to BLEND, but with slow pulses (i.e. for Polypectomies).
<b>MICRO</b>	Delicate Coagulation, with low sparking effect.
<b>M PULSED</b>	Delicate Micro Coagulation, with single pulses, which are adjustable from 0.1sec to 1sec. It is indicated for all the micro coagulations.
<b>FULG</b>	Macro Coagulation "Fulguration" with strong sparking effect. It is indicated to coagulate all the tissues, even under liquid, and also to perform very coagulating cuts.
<b>FULG PULSED</b>	Macro Coagulation "Fulguration", with fast pulses. Similar to FULG, but more delicate. It is indicated to reduce the smokes.

### BIPOLAR CURRENTS


<b>BIPOLAR</b>	Bipolar Coagulation, to be used with forceps, scissors, double needle electrodes and laparoscopic instruments.
----------------	--

## THE CURRENTS OF THE MODEL FOR BIPOLAR USE ONLY

<b>PURE</b>	Cut for the use in Open Surgery or Laparoscopy.
<b>P PULSED</b>	Pulsed Fast Cut, suitable to obtain very fine results, a minimum thermal effect and a reduction of the smokes.
<b>BLEND</b>	Coagulating Cut, with a higher thermal effect.
<b>MACRO</b>	Coagulation, to be used with forceps, scissors, double needle electrodes and laparoscopic instruments.
<b>M PULSED</b>	Pulsed Coagulation. Similar to MACRO, but more delicate and useful to reduce the smokes.
<b>MICRO</b>	Micro Coagulation, to be used with forceps, scissors, double needle electrodes and laparoscopic instruments.

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
For manual activation with model SU 140-MPC only, see: MPE/CMS, MP/CM

B730	ALSATOM SU 140/D-MPC, without accessories
B730/A	SET OF STANDARD ACCESSORIES FOR ALSATOM SU 140/D-MPC identical to B700/A, but with dual switch pedal D-STOP/P
B730/B	ditto, but with flexible, conductive silicone neutral plate NP/GP
B730/D	ditto, but for dental use, without EIP/9 and SEL/VI, with: EIP/S - hand cylinder neutral plate with cable 2.5 m and SEL/D - set of 8 dental electrodes

For manual activation see: MPE/CMS, MP/CM

B740	ALSATOM SU 140/BD-MPC, with: dual switch pedal D-STOP/P
------	--

## GENERAL LIST OF ACCESSORIES

### HANDLES

MPE/F	Autoclavable handle with cable mt. 2,5
MPE/CMS	Hand-switch handle (CUT/COAG), 100 times autoclavable with cable mt.3
MP/CM	Ditto, disposable (with knife or needle or ball electrode)

### Reusable NEUTRAL ELECTRODES

EIP/9	Stainless steel neutral electrode with cable mt.2,5
NP/GP	Ditto, flexible conductive silicone
EIP/S	Ditto, hand cylinder

### Disposable NEUTRAL ELECTRODES

CMS/E	Reusable connecting cable, mt. 2,5
EIP/SA	Adhesive dual electrode for adults (box 25pcs) (*)
EIP/SP	Ditto, paediatric (*)

### TROLLEYS AND FOOTSWITCHES

H10/AB	Two-shelves trolley
STOP/PN	Pneumatic, waterproof, explosion-proof single footswitch
D-STOP/P	Ditto, double

### SMOKE SUCTION UNITS

CS200	Smoke suction unit, with filter
-------	---------------------------------

### HANDLES AND NEEDLES FOR MICRO SURGERY

MLD/F	Autoclavable light handle, cable mt. 2
MID	Needle adaptor (for handles MPE series)
AID	Insulated needle - Ø 0.10 mm for depilation
SAD	Set of long and non insulated needles Ø 0.10mm for fine cut
SAD/1	Ditto, Ø 0.15 mm for fine cut, telangiectasia
SAD/2	Ditto, Ø 0.20 mm for fine cut, face depilation, telangiectasia of legs
SAD/3	Ditto, Ø 0.40 mm, for fine cut, microcoagulation

ACCESSORIES AND CABLES FOR BICOAGULATION, MICROBICOAGULATION, BICOAGULATION OF THE TURBINAL OR OF LARYNX, see Alsa catalogue: "Accessories for bicoagulation"

ELECTRODES FOR ELECTROSURGERY, FOR GYNAECOLOGY AND CONIZATION reusable or disposable, FOR DENTAL USE, see Alsa catalogue: "Monopolar electrodes"

ADAPTORS FOR MONOPOLAR CONNECTING CABLES (e.g. accessories for Endoscopy) OF BIPOLAR ACCESSORIES (made by other manufacturers), see Alsa catalogue: "Adapters"

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